

## Product datasheet

# Recombinant Hepatitis C Virus Hepatitis C Virus NS4 protein ab49026

### 製品の詳細

---

製品名	Recombinant Hepatitis C Virus Hepatitis C Virus NS4 protein
生理活性	Strongly reacts with human HCV positive serum.
精製度	> 95 % SDS-PAGE. Purity of proteins is evaluated by SDS-PAGE
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Hepatitis C virus
配列	MNRLI AFASRGNHVS PTHYVPESDA AARVTQI
領域	1916 to 1947

### 特性

---

Our **Abpromise guarantee** covers the use of **ab49026** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA Flow Cytometry SDS-PAGE
----------	---

製品の状態 Liquid

備考 Strongly reacts with human HCV positive serum.

Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.

It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

### 前処理および保存

---

## 保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

pH: 7.2

Constituents: 0.2% Triton-X-100, 9% Urea, 0.395% Tris HCl, 50% Glycerol (glycerin, glycerine)

## 関連情報

---

### 関連性

NS4 A/B are two of the seven nonstructural (NS) proteins making up the HCV polyprotein. The NS proteins recruit the viral genome into an RNA replication complex, which is associated with rearranged cytoplasmic membranes. NS4A acts as a cofactor with the NS3 serine protease and stabilizes its folding. The NS3-NS4A complex is essential for the activation of the latter and allows membrane anchorage of NS3. HCV is a positive, single-stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of about 3,000 amino acids. HCV is responsible for a large proportion of worldwide chronic viral hepatitis. Most of these infections develop into chronic hepatitis, which often progresses to liver cirrhosis and hepatocellular carcinoma. At present, (unlike hepatitis A and B), there is no vaccine to prevent hepatitis C infection. The hepatitis C virus (HCV) nonstructural protein 4B (NS4B) is a relatively hydrophobic 27-kDa protein. The 4A protein has a molecular weight of 6 kDa.

### 細胞内局在

Endoplasmic reticulum; multi-pass membrane protein.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors