

Recombinant HIV1 + HIV2 envelope antigen protein ab48966

製品の詳細

製品名	Recombinant HIV1 + HIV2 envelope antigen protein
生理活性	Detects all HIV1,2 infected individuals responding to HIV1,2 envelope proteins.
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant

特性

Our **Abpromise guarantee** covers the use of **ab48966** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot
製品の状態	Liquid
備考	Detects all HIV1,2 infected individuals responding to HIV1,2 envelope proteins.

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 9% Urea, 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol (glycerin, glycerine)
------------	--

関連情報

関連性	Inhibitors targeting HIV1 viral envelope proteins are used as antiretroviral drugs. The surface protein gp120 attaches the virus to the host lymphoid cell by binding to the primary receptor CD4. This interaction induces a structural rearrangement creating a high affinity binding site for a chemokine coreceptor like CXCR4 and/or CCR5. The transmembrane protein gp41 acts as a class I viral fusion protein.
-----	--

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors