

### Recombinant HIV1 p17 protein ab63981

画像数 1

#### 製品の詳細

製品名	Recombinant HIV1 p17 protein
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant

#### 特性

Our **Abpromise guarantee** covers the use of **ab63981** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	ab63981 has been derived from infectious plasmid clone pNL4-3 and belongs to Group B of HIV-1.

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. pH: 6 Constituents: 0.078% Beta mercaptoethanol, 0.316% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.29% Sodium chloride
------------	--

#### 関連情報

関連性	HIV1 p17 is the matrix protein of the Gag polyprotein which performs highly complex orchestrated tasks during the assembly, budding, maturation, and infection stages of the viral replication cycle. During viral assembly, the proteins form membrane associations and self-associations that ultimately result in budding of an immature virion from the infected cell. Gag precursors also function during viral assembly to selectively bind and package two plus strands of genomic RNA.
-----	--

## 画像



Polyacrylamide gel electrophoresis of ab63981.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors