

Recombinant Human ASC1 protein ab152970

画像数 1

製品の詳細

製品名 Recombinant Human ASC1 protein

発現系 Wheat germ

タンパク質長 Full length protein

Animal free No

由来 Recombinant

生物種 Human

配列

```
MAVAGAVSGEPLVHWCTQQLRKTFGLDVSEEEIIQYVLSIESA
EEIREYVT
DLLQGNKGGKGFIEELITKWQKNDQELISDPLQQCFKKDEI
LDGQKSGD
HLKRGRKKGRNRQEVPAFTEPDTTAEVKTPFDLAKAQENSNS
VKKKTKFV
NLYTREGQDR LAVLLPGRHPCDCLGQKHKLINNCLICGRIVC
EQEGSGPC
LFCGTLVCTHEEQDILQRDSNKSQKLLKLLMSGVENSNGKVDI
STKDLLPH
QELRIKSGLEKAIKHKDKLLEFDRTSIRRTQVIDDESDFAS
DSNQWLSK
LERETLQKREELREL RHASRLSKKVTIDFAGRKILEEENSL
AEYHSRLD
ETIQAIANGTLNQPLTKLDRSSEELGVLVNPVQSPQWV
DHTGAASQ
KKAFRSSGFGLEFNSFQHLRIQDQEFQEGFDGGWCLSVHQP
WASLLVRG
IKRVEGRSWYTPHRGRLWIAATAKKPSPQEVSELQATYRLLR
GKDVEFPN
DYPSGCLLGCVLDLIDCLSQKQFKEQFPDISQESDSPFVFICK
NPQEMVVK FPIKGNPKIWKLDSKIHQGAKKGLMKQNKAV
```

領域 1 to 581

タグ GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab152970** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Western blot

ELISA

製品の状態 Liquid

備考

前処理および保存

保存方法および安定性 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00


Constituents: 0.31% Glutathione, 0.79% Tris HCl

関連情報

機能 Transcription coactivator of nuclear receptors which functions in conjunction with CBP-p300 and SRC-1 and may play an important role in establishing distinct coactivator complexes under different cellular conditions. Plays a pivotal role in the transactivation of NF-kappa-B, SRF and AP1. Acts as a mediator of transrepression between nuclear receptor and either AP1 or NF-kappa-B. Plays a role in androgen receptor transactivation and in testicular function.

細胞内局在 Nucleus. Cytoplasm. Cytoplasmic under conditions of serum deprivation.

画像



ab152970 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human ASC1 protein
(ab152970)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors