

### Recombinant Human BCMA protein (Tagged) (Biotin) ab269984

#### 画像数 1

#### 製品の詳細

---

製品名	Recombinant Human BCMA protein (Tagged) (Biotin)
精製度	> 73 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	<b>Q02223</b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MLQMAGQCSQNEYFDSLLHACIPCQLRCSSTPPLTCQRYCN ASVTNSVK GTNA
予測される分子量	34 kDa including tags
実際の分子量	37 kDa including tags
分子量情報	34 kDa + Glycans
領域	1 to 54
タグ	Avi tag C-Terminus , Fc tag C-Terminus
配列の追加情報	a.a. 1-54 fused to Fc region of Human IgG1. NM_001192
標識	Biotin

#### 特性

---

Our **Abpromise guarantee** covers the use of **ab269984** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

#### 前処理および保存

---

保存方法および安定性 Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40


Constituents: 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine), 0.13% Sodium phosphate

Protein may be diluted to 100 µg/ml in PBS + glycerol and stored at -80°C.

## 関連情報

<b>機能</b>	Receptor for TNFSF13B/BLYS/BAFF and TNFSF13/APRIL. Promotes B-cell survival and plays a role in the regulation of humoral immunity. Activates NF-kappa-B and JNK.
<b>組織特異性</b>	Expressed in mature B-cells, but not in T-cells or monocytes.
<b>関連疾患</b>	Note=A chromosomal aberration involving TNFRSF17 is found in a form of T-cell acute lymphoblastic leukemia (T-ALL). Translocation t(4;16)(q26;p13) with IL2.
<b>配列類似性</b>	Contains 1 TNFR-Cys repeat.
<b>細胞内局在</b>	Cell membrane. Endomembrane system. Perinuclear Golgi-like structures.

## 画像



SDS-PAGE analysis of 4 µg ab269984 (Lane 1). 4-20% gel with Coomassie staining.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors