

Product datasheet

Recombinant Human BTN3A1 protein (Tagged) ab271422

画像数 2

製品の詳細

製品名	Recombinant Human BTN3A1 protein (Tagged)
精製度	>= 90 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	<u>O00481</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	FSVLGPSGPI LAMVGEDADL PCHLFPTMSA ETMELKQVSS SLRQVVNVYA DGKEVEDRQS APYRGRSIL RDGITAGKAA LRIHNVTASD SGKYLCYFQD GDFYEKALVE LKVAALGSDL HVDVKGYKDG GIHLECRSTG WYPQPQIQWS NNKGENIPTV EAPVVADGVG LYAVAASVIM RGSSGEGVSC TIRSSLLGLE KTASISIADP FFRSAQRWIA ALAG
予測される分子量	53 kDa including tags
領域	31 to 254
タグ	Avi tag C-Terminus , Fc tag C-Terminus
配列の追加情報	Fc region of Human IgG1. Genbank: NM_007048

特性

Our **Abpromise guarantee** covers the use of **ab271422** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

関連情報

組織特異性

Highly expressed in heart, pancreas and lung, Moderately expressed in placenta, liver and muscle.

配列類似性

Belongs to the immunoglobulin superfamily. BTN/MOG family.


Contains 1 B30.2/SPRY domain.

Contains 2 Ig-like V-type (immunoglobulin-like) domains.

細胞内局在


Membrane.

画像



SDS-PAGE analysis of ab271422.

SDS-PAGE - Recombinant Human BTN3A1 protein
(Tagged) (ab271422)



All lanes :

Lane 1 : ab271422

Lane 2 : MW markers

Predicted band size: 58 kDa

Western blot - Recombinant Human BTN3A1 protein
(Tagged) (ab271422)

Western blot analysis of ab271422.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors