

# Recombinant Human Cbl-c protein (Tagged) ab271396

画像数 1

### 製品の詳細

製品名	Recombinant Human Cbl-c protein (Tagged)
精製度	>= 90 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	<b><u>Q9UIV8</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	QWEEARALGR AVRMLQRLEE QCVDPRLSVS PPSLRDLLPR TAQLLREVAH SRRAAGGGGP GGPGGSGDFL LIYLANLEAK SRQVAALLPP RGRRSANDEL FRAGSRLRRQ LAKLAIIFSH MHAELHALFP GGKYCGHMYQ LTKAPAHTFW RESCGARCVL PWAEFESLLG TCHPVEPGCT ALALRTTIDL TCSGHVSIFE FDVFTRLFQP WPTLLKNWQL LAVNHPGYMA FLTYDEVQER LQACRDKPGS YIFRPSCTRL GQWAIGYVSS DGSILQTIPA NKPLSQVLE GQKDGFFLYP DGKTHNPDLT ELGQAEPQQR IHVSEEQLQL YWAMDSTFEL CKICAESNKD VKIEPCGHLL CSCCLAAWQH SDSQTCPFQR CEIKGWEAVS IYQF
予測される分子量	45 kDa including tags
領域	11 to 404
タグ	DDDDK tag N-Terminus
配列の追加情報	Genbank: NM_012116

### 特性

Our **Abpromise guarantee** covers the use of **ab271396** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

## 製品の状態

Liquid

## 前処理および保存

### 保存方法および安定性

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 0.04% Tween, 20% Glycerol (glycerin, glycerine)

100 µg/mL DDDDK peptide.

## 関連情報


### 関連性

Cblc belongs to a family of ubiquitin protein ligases (E3s) that negatively regulate signaling by targeting activated tyrosine kinases for degradation. Cblc is the most recently cloned member of the Cbl proteins and is expressed only in epithelial cells (the other Cbl proteins are ubiquitous). Cblc, like the other mammalian Cbl proteins, can ubiquitinate the activated EGFR and target it for degradation. Cblc knock out mice show no obvious phenotype. Thus, the physiological role of Cblc is not known.

### 細胞内局在

Nuclear

## 画像



SDS-PAGE analysis of ab271396.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors