

Recombinant Human CD164 protein (His tag) ab227398

画像数 1

製品の詳細

製品名	Recombinant Human CD164 protein (His tag)
精製度	> 85 % SDS-PAGE. Affinity purified
エンドキシン・レベル	< 1.000 Eu/μg
発現系	Baculovirus infected insect cells
アクセッション番号	Q04900
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	ADPDKNTTQH PNVTTLAPIS NVTSAVPTSL PLVTTTPAPET CEGRNSCVSC FNVSVVNTTC FWIECKDESY CSHNSTVSDC QVGNTTDFCS VSTATPVPTA NSTAKPTVQP SPSTTSKTVT TSGTTNNTVT PTSQPVRKST FDHHHHHH
予測される分子量	16 kDa including tags
領域	24 to 162
タグ	His tag C-Terminus
配列の追加情報	Extracellular domain. NP_006007

特性

Our **Abpromise guarantee** covers the use of **ab227398** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------	--

pH: 7.40

Constituents: PBS, 10% Glycerol (glycerin, glycerine)

関連情報

関連性 CD164 is a member of the CD164 family of proteins and is highly N- and O-glycosylated. It has a Ser-Gly motif that may serve as the site of attachment of a glycosaminoglycan side chain. It contains sialic acid. CD164 may be associated with carcinoma and is probably a mucin.

細胞内局在 Secreted

画像



15% SDS-PAGE - 3 µg Recombinant Human CD164 protein (His tag) (ab227398).

Migrates as 28-57KDa (SDS-PAGE under reducing conditions).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors