

Recombinant Human CD1a protein ab114359

画像数 1

製品の詳細

製品名	Recombinant Human CD1a protein
発現系	Wheat germ
アクセッション番号	P06126
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MWNWLKEPLSFHVIWIASFYNSHWKQNLVSGWLSLQHTHTWD SNSSTIVF LWPWSRGNFSNEEWKELETFRIRTIKRSFEGIRRYAHELQFE YPFEIQVT GGCELHSGKVSFSLQLAYQGSDFVSFQNSWLPPVAGNMA KHFCKVLN QNHENDITHNLLSDTCPRFILGLLDAGKAHLQRQVKPEAWL SHGPPSPG GHLQLVCHVSGFYKPVWMMRGEQEQGTQRGDILPSADG TWSLRATL EVAAGEAADLSCRVKHSSLEGQDIVLYWEHSSVGFIIILAVI VPLLLLIG LALWFRKRCFC
予測される分子量	60 kDa including tags
領域	1 to 312

特性

Our **Abpromise guarantee** covers the use of **ab114359** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid

保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能

Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells.

組織特異性

Expressed on cortical thymocytes, epidermal Langerhans cells, dendritic cells, on certain T-cell leukemias, and in various other tissues.


配列類似性

Contains 1 Ig-like (immunoglobulin-like) domain.

細胞内局在

Cell membrane. Endosome membrane. Subject to intracellular trafficking between the cell membrane and endosomes. Localizes to cell surface lipid rafts.

画像



12.5% SDS-PAGE showing ab114359 at approximately 59.95kDa stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors