

Product datasheet

Recombinant Human CD58 protein (denatured) ab187438

画像数 1

製品の詳細

製品名	Recombinant Human CD58 protein (denatured)
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	<u>P19256</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMGFSQIYGVVYGNVTFHVPS NVPLKEVL WKKQKDKVAELENSEFRAFSSFKNRVYLDTVSGSLTIYNLTS SDEDEYEM ESPNITDTMKFFLYVLES LPSPTLTCAL TNGSIEVQCM IPEH YNSHRGLI MYSWDCPMEQCKRNSTSIYFKMENDLPQKIQCTLSNPLFNTT SSIILTTC IPSSGHSRHR
予測される分子量	24 kDa including tags
領域	29 to 215
タグ	His tag N-Terminus
配列の追加情報	Extracellular domain (NP_001770).

特性

Our **Abpromise guarantee** covers the use of **ab187438** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	Previously labelled as LFA3.

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.4

Constituents: 10% Glycerol (glycerin, glycerine), 90% Phosphate Buffer

関連情報

機能

Ligand of the T-lymphocyte CD2 glycoprotein. This interaction is important in mediating thymocyte interactions with thymic epithelial cells, antigen-independent and -dependent interactions of T-lymphocytes with target cells and antigen-presenting cells and the T-lymphocyte rosetting with erythrocytes. In addition, the LFA-3/CD2 interaction may prime response by both the CD2+ and LFA-3+ cells.


配列類似性

Contains 1 Ig-like C2-type (immunoglobulin-like) domain.

細胞内局在

Cell membrane.

画像



15% SDS-PAGE analysis of 3 µg of ab187438.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors