

Recombinant Human CD81 protein ab152267

画像数 1

製品の詳細

製品名	Recombinant Human CD81 protein
発現系	Wheat germ
アクセッション番号	P60033
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	GGVILGVALWLRHDPQTTNLLYLELGDKPAPNTFYVGIYILI AVGAVMMF VGFLGCGAIQESQCLLGTFFTCLVILFACEVAAGIWGFVNK DQIAKDVK QFY
予測される分子量	37 kDa including tags
領域	25 to 127
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab152267** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot

製品の状態 Liquid

備考

前処理および保存


保存方法および安定性 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

関連情報

機能	May play an important role in the regulation of lymphoma cell growth. Interacts with a 16-kDa Leu-13 protein to form a complex possibly involved in signal transduction. May acts as the viral receptor for HCV.
組織特異性	Hematolymphoid, neuroectodermal and mesenchymal tumor cell lines.
関連疾患	Defects in CD81 are the cause of immunodeficiency common variable type 6 (CVID6) [MIM:613496]; also called antibody deficiency due to CD81 defect. CVID6 is a primary immunodeficiency characterized by antibody deficiency, hypogammaglobulinemia, recurrent bacterial infections and an inability to mount an antibody response to antigen. The defect results from a failure of B-cell differentiation and impaired secretion of immunoglobulins; the numbers of circulating B cells is usually in the normal range, but can be low.
配列類似性	Belongs to the tetraspanin (TM4SF) family.
翻訳後修飾	Not glycosylated.
細胞内局在	Membrane.

画像



ab152267 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors