

### Recombinant Human CD89 protein ab151878

#### 製品の詳細

製品名	Recombinant Human CD89 protein
精製度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE. ab151878 was lyophilised from a 0.2 µM filtered solution.
エンドトキシン・レベル	< 1.000 Eu/µg
発現系	HEK 293 cells
アクセッション番号	<b><u>P24071</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	QEGDFMPFISAKSSPVIPLDGSVKIQCAIREAYLTQLMII KNSTYREI GRRLKFWNETDPEFVIDHMDANKAGRYQCQYRIGHYRFRYSD TLELVVTG LYGKPFLSADRGLVLMPGENISLTCSSAHIPFDRFSLAKEGE LSLPQHQS GEHPANFSLGPVDLNVSGIYRCYGWYNRSPYLWSFPSNALEL VVTDSIHQ DYTTQNVDDHHHHHH
予測される分子量	25 kDa including tags
領域	22 to 227
タグ	His tag C-Terminus

#### 特性

Our **Abpromise guarantee** covers the use of **ab151878** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	HPLC

製品の状態	Lyophilized
-------	-------------

#### 前処理および保存

<b>保存方法および安定性</b>	Shipped at 4°C. Store at -20°C or -80°C. pH: 7.20 Constituents: 99% Phosphate Buffer, 0.88% Sodium chloride
<b>再構成</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilised protein in 1X PBS. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. For long term storage aliquot and store at < -20°C.
<b>関連情報</b>	
<b>機能</b>	Binds to the Fc region of immunoglobulins alpha. Mediates several functions including cytokine production.
<b>組織特異性</b>	Isoform A.1, isoform A.2 and isoform A.3 are differentially expressed between blood and mucosal myeloid cells. Isoform A.1, isoform A.2 and isoform A.3 are expressed in monocytes. Isoform A.1 and isoform A.2 are expressed in alveolar macrophages; however only one isoform is expressed at alveolar macrophages surfaces.
<b>配列類似性</b>	Contains 2 Ig-like C2-type (immunoglobulin-like) domains.
<b>細胞内局在</b>	Secreted and Cell membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors