

Recombinant Human CDKN2A/p16INK4a protein ab84075

2 References 画像数 1

製品の詳細

製品名	Recombinant Human CDKN2A/p16INK4a protein
精製度	> 95 % SDS-PAGE. Purity is greater than 95% by SDS-PAGE gel and HPLC analyses.
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MEPAAGSSME PSADWLATAA ARGRVEEVRA LLEAGALPNA PNSYGRRPIQ VMMMG SARVA ELLLLHGAEP NCADPATLTR PVHDAAREGF LDTLVVLHRA GARLDVRDAW GRLPVDLAEE LGHRDVARYL RAAAGGTRGS NHARIDAAEG PSDIPD
予測される分子量	17 kDa
領域	1 to 156

特性

Our **Abpromise guarantee** covers the use of **ab84075** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	Endotoxin level is less than 0.1 ng per μg (1EU/μg).

前処理および保存


保存方法および安定性	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. Constituents: 0.00268% PBS, 0.0154% DTT
再構成	Reconstitute with 1x PBS, pH 8.0. to a concentration of 0.1-1.0 mg/ml.

関連情報

細胞内局在

Cytoplasmic and Nuclear

画像



Standard Curve for CDKN2A/p16INK4a (Analyte:

CDKN2A/p16INK4a protein (Human) (ab84075)); dilution range

1pg/ml to 1µg/ml using Capture Antibody **Mouse monoclonal**

[DCS50.1] to CDKN2A/p16INK4a (ab16123) at 1µg/ml and

Detector Antibody **Rabbit polyclonal to CDKN2A/p16INK4a**

(ab7962) at 0.5µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors