

Recombinant Human CEAP protein ab126912

画像数 1

製品の詳細

製品名	Recombinant Human CEAP protein
精製度	> 85 % SDS-PAGE. ab126912 is purified by conventional chromatography techniques to > 85% by SDS - PAGE.
発現系	Escherichia coli
アクセッション番号	<u>Q6QNY1</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHH SSGLVPRGSH MGSMAAAAE GVLATRSDEP ARDDAAVETA EEAKEPAEAD ITELCRDMFS KMATYLTGEL TATSEYKLLNMNKLTSLK YLEMKDIAIN ISRNLDLNQ KYAGLQPYLD QINVIEEQVA ALEQAAYKLD AYSKKLEAKY KKLEKR
予測される分子量	19 kDa including tags
領域	1 to 142
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab126912** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as BLOC1S2

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

関連情報

機能

May play a role in cell proliferation. The BLOC-1 complex is required for normal biogenesis of lysosome-related organelles, such as platelet dense granules and melanosomes. Plays a role in intracellular vesicle trafficking.

組織特異性

Isoform 1 and isoform 2 are widely expressed. Expressed in various malignant tumor tissues (at protein level).


配列類似性

Belongs to the BLOC1S2 family.

細胞内局在

Cytoplasm > cytoskeleton > centrosome. Localizes to the centrosomes in a microtubule-dependent manner.

画像



15% SDS-PAGE showing ab126912 (3 µg). Please note: Molecular weight on SDS-PAGE will appear higher.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors