

Recombinant Human Chromogranin B protein ab93930

画像数 1

製品の詳細

製品名	Recombinant Human Chromogranin B protein
精製度	> 80 % SDS-PAGE. Purified using conventional chromatography techniques
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
領域	21 to 677

特性

Our **Abpromise guarantee** covers the use of **ab93930** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
------------	---

関連情報

機能	Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active peptides.
組織特異性	Expressed in the adrenal medulla, and in pheochromocytoma. Not expressed in liver.
配列類似性	Belongs to the chromogranin/secretogranin protein family.


翻訳後修飾

Extensively processed by limited proteolysis at conserved basic residues. Alternative processing are seen in different tissues.

細胞内局在

Secreted. Neuroendocrine and endocrine secretory granules.

画像



15% SDS-PAGE showing ab93930 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors