

Recombinant Human coronavirus MERS Spike glycoprotein (S1) (His tag) ab270821

画像数 1

製品の詳細

製品名	Recombinant Human coronavirus MERS Spike glycoprotein (S1) (His tag)
精製度	> 95 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human coronavirus
領域	56 to 295
タグ	His tag C-Terminus
配列の追加情報	6xHis tag
製品の詳細	Recombinant Human coronavirus MERS Spike glycoprotein (His tag)

特性


Our **Abpromise guarantee** covers the use of **ab270821** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Immunoassay SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped on Dry Ice. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7 Preservative: 0.05% Sodium azide Constituents: 0.44% L-Arginine, PBS
------------	---



SDS-PAGE analysis of ab270821.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors