

# Recombinant Human coronavirus NL63 Nucleoprotein ab270843

画像数 1

### 製品の詳細

製品名	Recombinant Human coronavirus NL63 Nucleoprotein
精製度	> 95 % SDS-PAGE. Purified by affinity chromatography (nickel column).
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human coronavirus
タグ	His tag C-Terminus
配列の追加情報	Human Coronavirus NL63 Nucleoprotein C-terminal region (130 amino acids).

### 特性

Our **Abpromise guarantee** covers the use of **ab270843** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

### 前処理および保存

保存方法および安定性	Shipped on Dry Ice. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7 Preservative: 0.05% Sodium azide Constituent: PBS
------------	--

### 画像



Coomassie-stained SDS-PAGE showing ab270843.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors