

## Product datasheet

# Recombinant Human COUP TF1 protein ab152739

画像数 1

### 製品の詳細

製品名	Recombinant Human COUP TF1 protein
発現系	Wheat germ
アクセッション番号	<b>P10589</b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MAMVVSSWRDPQDDVAGGNPGGPNPAAQAARGGGGGAGEQQQ QAGSGAPH TPQTPGQPGAPATPGTAGDKGQGGPPGSGQSQQHIECVVCGDK SSGKHYGQ FTCEGCKSFFKRSVRRNLTYTCRANRNCPIDQHHRNQCQYCR LKKCLKVG MRREAVQRGRMPPTQPNPGQYALTNGDPLNGHCYLSGYISLL LRAEPYPT SRYGSQCMQPNNIMGIENICELAARLLFSAVEWARNIPFFPD LQITDQVS LLRLTWSELFVLNAAQCSMPLHVAPLLAAAGLHASPMSADRV VAFMDHIR IFQEQVEKALKALHVDSAEYSCLKAIVLFTSDACGLSDAAHIE SLQEKSQC ALEEYVRSQYPNQPSRFGKLLLRLPSLRTVSSSVIEQLFFVR LVGKTPIE TLIRDMLLSGSSFNWPYMSIQCS
予測される分子量	73 kDa including tags
領域	1 to 423

### 特性

Our **Abpromise guarantee** covers the use of **ab152739** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE

	ELISA
製品の状態	Liquid
備考	


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------	--

#### 関連情報

機能	Coup (chicken ovalbumin upstream promoter) transcription factor binds to the ovalbumin promoter and, in conjunction with another protein (S300-II) stimulates initiation of transcription. Binds to both direct repeats and palindromes of the 5'-AGGTCA-3' motif. Represses transcriptional activity of LHCG.
配列類似性	Belongs to the nuclear hormone receptor family. NR2 subfamily. Contains 1 nuclear receptor DNA-binding domain.
細胞内局在	Nucleus.

#### 画像



12.5% SDS-PAGE analysis of ab152739 stained with Coomassie Blue.

SDS-PAGE - Recombinant Human COUP TF1 protein (ab152739)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors