

Recombinant Human CrkL protein ab116208

画像数 1

製品の詳細

製品名	Recombinant Human CrkL protein
精製度	> 90 % SDS-PAGE. ab116208 was purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>P46109</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMSSARFDSSDRSAWYMGPVSRQEAQTRLQGQRHGMFLVRDSSTCPGDYVLSVSENSRVSHYIINSLPNRRFKIGDQEFDLPALLEFYKIHLYLDTTTLIEPAPRYPSPPMGSVSAPNLPTAEDNLEYVRTLYDFPGNDAEDLPFKKGEILVIEKPEEQWWSARNKDGRVGMIPVPYVEKLVRSSPHGKHGNRNSNSYGIPEPAHAYAQPQTTTLPVAVSGSPGAAITPLPSTQNGPVFAKAIQKRVPAYDKTALALEVGDIVKVTMNIINQWEGEVN GRKGLFPFTHVKIFDPQNPENE
予測される分子量	36 kDa including tags
領域	1 to 303
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab116208** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF

製品の状態

Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能

May mediate the transduction of intracellular signals.


配列類似性

Belongs to the CRK family.

Contains 1 SH2 domain.

Contains 2 SH3 domains.

画像



15% SDS-PAGE analysis of 3 µg ab116208.

SDS-PAGE - Recombinant Human CrkL protein
(ab116208)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors