

### Recombinant Human E2F8 protein ab153513

画像数 1

#### 製品の詳細

製品名	Recombinant Human E2F8 protein
発現系	Wheat germ
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MAQLAAICKMQLEEQSSESQRQVKVQLARSGPCKPVAPLDPP VNAEMELT APSLIQPLGMVPLIPSPSSAVPLILPQAPSGPSYAIYLQPT QAHQSVTP PQGLSPTVCTTHSSKATGSKDSTDATTEKAANDTSKASASTR PGSLLPAP ERQGAKSRTREPAGERGSKRASMLEDSGSKKKFKEDLKLEN VSATLFPS GYLIPLTQCSSLGAESILSGKENSALSHPNHRIYSSPIAGVI PVTSSSELT AVNFPSFHVTPLKLMVSPTSVAAPVGNPALASSHPVPIQN PSSAIVNF TLQHLGLISPNVQLSASPGSGIVPVSPRIESVNVAPENAGTQ QGRATNYD SPVPGQSQPNGQSVAVTGAQQPVPVTPKGSQLVAESFFRTPG GPTKPTSS SCMDFEGANKTSLGTLFVPQRKLEVSTEDVH
領域	1 to 431
タグ	GST tag N-Terminus

#### 特性

Our **Abpromise guarantee** covers the use of **ab153513** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	ELISA
製品の状態	Liquid

## 備考

### 前処理および保存

#### 保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

### 関連情報

#### 機能

Along with E2F7, inhibitor of E2F-dependent transcription. Binds DNA independently of DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3'. Directly represses E2F1 transcription. Appears to regulate a subset of E2F-dependent genes whose products are required for normal cell cycle progression.

#### 配列類似性

Belongs to the E2F/DP family.


#### ドメイン

Both DNA-binding domains are required for DNA-binding and are proposed to form an intramolecular structure that is similar to the winged helix structure of the E2F-DP heterodimer.

#### 細胞内局在

Nucleus.

### 画像



ab153513 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human E2F8 protein

(ab153513)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors