

Recombinant Human EMR1/ADGRE1 (F4/80) protein ab131740

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant Human EMR1/ADGRE1 (F4/80) protein |
| 発現系 | Wheat germ |
| アクセッション番号 | <u>Q14246-2</u> |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MRGFNLLLFWGCCVMHSWEGHIRPTRKPNKGNCRDSTLCP AYATCTNT VDSYYCACKQGFLSSNGQNHFKDPGVRCKDIDECSSQPQPCG PNSSCKNL SGRYKCSCLDGFSSPTGNDWVPGKPGNFSCDINECLTSSVC PEHSDCVN SMGSYSCSCQVGFISRNSTCEDVDECADPRACPEHATCNNTV GNYSCFCN PGFESSGHL SFQGLKASCEDIDECTEMCPINSTCTNTPGSY FCTCHPGF APSNGQLNFTDQGV ECRDIDECRQDPSTCGPNSICTNALGSY SCGCIAGF HPNPEGSQKDNFSCQRVLFKCKEDVIPDNKQIQQCQEGTAV KPAYVSFC AQINNIFSVLDKVCENKTTVVSLKNTTESFVPVLKQISTWTK FTKEETSS LATVFLESVESMTLASFVKPSANITPAVRTEYLDIESKVINK ECSEENV LDLVAKGDKMKIGCSTIEESESTETTGVAFVSFVGMESVLNE RFFKDHQA PLTTSEIKLKMNSRVVGGIMTGEKKDGFSDPIIYTLENIQPK QKFERPIC VSWSTDVKGGRWTSFGCVILEASETYTICSCNQMANLAIIMA SGELTMGC AIIAGFLHYLFLACFFWMLVEAVILFLMVRNLKVVNYFSSRN IKMLHICA FGYGLPMLVVVISASVQPQGYGMHNRCWLNNTETGFIWSFLGP VCTVIVIN |

SLLLTWTLWILRQLSSVNAEVSTLKDTRLLTFKAFAQLFIL
GCSWVLGI
FQIGPVAGVMAYLFTIINSLQGAFIFLIHCLLNGQVREEYKR
WITGKTKP SSQSQTSRILLSSMPSASKTG

| | |
|----------|---------------------------------------|
| 予測される分子量 | 117 kDa including tags |
| 領域 | 1 to 821 |
| 製品の詳細 | Recombinant Human EMR1/ADGRE1 protein |

特性

Our **Abpromise guarantee** covers the use of **ab131740** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|--------------|
| アプリケーション | SDS-PAGE |
| | Western blot |
| | ELISA |

| | |
|-------|--------|
| 製品の状態 | Liquid |
|-------|--------|

備考

前処理および保存

| | |
|------------|--|
| 保存方法および安定性 | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl |
|------------|--|

関連情報

| | |
|-------|---|
| 機能 | Orphan receptor involved in cell adhesion and probably in cell-cell interactions specifically involving cells of the immune system. May play a role in regulatory T-cells (Treg) development. |
| 組織特異性 | Expression is restricted to eosinophils. |
| 配列類似性 | Belongs to the G-protein coupled receptor 2 family. Adhesion G-protein coupled receptor (ADGR) subfamily. Contains 6 EGF-like domains. Contains 1 GPS domain. |
| 細胞内局在 | Cell membrane. |

画像



12.5% SDS-PAGE analysis of ab131740, stained with Coomassie Blue.

SDS-PAGE - Recombinant Human EMR1/ADGRE1 protein (ab131740)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors