

Recombinant Human EPCR/CD201 protein ab151362

製品の詳細

製品名	Recombinant Human EPCR/CD201 protein
精製度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
エンドトキシン・レベル	< 1.000 Eu/μg
発現系	HEK 293 cells
アクセッション番号	<u>Q9UNN8</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	SQDASDGLQRLHMLQISYFRDPYHVWYQGNASLGGHLTHVLE GPDTNTTI IQLQPLQEPESWARTQSGLQSYLLQFHGLVRLVHQERTLAFP LTIRCFLG CELPPEGSRAHVFFEVAVNGSSFVSRPERALWQADTQVTSG VVTFTLQQ LNAYNRTRYELREFLEDTCVQYVQKHISAENTKGSQTSRSYT SVDHHHHH H
予測される分子量	23 kDa including tags
領域	18 to 210
タグ	His tag C-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab151362** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE HPLC
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -80°C. pH: 7.20 Constituents: 99% Phosphate Buffer, 0.88% Sodium chloride
再構成	Dissolve the lyophilized protein in 1X PBS. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
関連情報	
機能	Binds activated protein C. Enhances protein C activation by the thrombin-thrombomodulin complex; plays a role in the protein C pathway controlling blood coagulation.
組織特異性	Expressed strongly in the endothelial cells of arteries and veins in heart and lung, less intensely in capillaries in the lung and skin, and not at all in the endothelium of small vessels of the liver and kidney.
翻訳後修飾	N-glycosylated. A soluble form exists; probably released by a metalloprotease. Seems to have the same activity as the membrane-bound form.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors