

# Recombinant human FIST protein ab85601

画像数 2

### 製品の詳細

---

製品名	Recombinant human FIST protein
生理活性	The Specific activity of ab85601 was determined to be 48 nmol/min/mg.
精製度	> 85 % Densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>Q9H422-1</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	163 to 562

### 特性

---

Our **Abpromise guarantee** covers the use of **ab85601** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot Functional Studies
製品の状態	Liquid
備考	<b>ab64311</b> (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity Previously labelled as HIPK3.

### 前処理および保存

---


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.00174% PMSF, 0.82% Sodium phosphate, 0.00308% DTT, 25% Glycerol (glycerin, glycerine), 1.74% Sodium chloride
------------	--

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 関連情報

機能	Seems to negatively regulate apoptosis by promoting FADD phosphorylation. Enhances androgen receptor-mediated transcription. May act as a transcriptional corepressor for NK homeodomain transcription factors.
組織特異性	Overexpressed in multidrug resistant cells. Highly expressed in heart and skeletal muscle, and at lower levels in placenta, pancreas, brain, spleen, prostate, thymus, testis, small intestine, colon and leukocytes. Not found in liver and lung.
配列類似性	Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family. HIPK subfamily. Contains 1 protein kinase domain.
翻訳後修飾	Autophosphorylated, but autophosphorylation is not required for catalytic activity. May be sumoylated.
細胞内局在	Cytoplasm. Nucleus.

## 画像



Functional Studies - Recombinant human FIST protein (ab85601)

The Specific activity of ab85601 was determined to be 48 nmol/min/mg.



SDS-PAGE - Recombinant human FIST protein (ab85601)

SDS-PAGE showing ab85601 at approximately 49kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors