

Recombinant Human GTR protein ab208317

画像数 1

製品の詳細

製品名	Recombinant Human GTR protein
精製度	> 95 % SDS-PAGE. Affinity purified
エンドキシン・レベル	< 1.000 Eu/μg
発現系	Insect cells
アクセッション番号	Q9Y5U5
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	QRPTGGPGCGPGRLLLLGTGTDARCCRVHTTRCCRDYPGEECC SEWDCMCV QPEFHCGDPCCTTCRHHPCPPGQGVQSQGKFSFGFQCIDCAS GTFSGGHE GHCKPWT DCTQFGFLT VFPGNKTHNAVCPGSPPAEPLHHH HHH
予測される分子量	16 kDa including tags
領域	26 to 162
タグ	His tag C-Terminus
配列の追加情報	NP_004186.

特性

Our **Abpromise guarantee** covers the use of **ab208317** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as TNFRSF18

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 10% Glycerol (glycerin, glycerine), 90% PBS

関連情報

機能

Receptor for TNFSF18. Seems to be involved in interactions between activated T-lymphocytes and endothelial cells and in the regulation of T-cell receptor-mediated cell death. Mediated NF-kappa-B activation via the TRAF2/NIK pathway.

組織特異性

Expressed in lymph node, peripheral blood leukocytes and weakly in spleen.


配列類似性

Contains 3 TNFR-Cys repeats.

細胞内局在

Secreted and Cell membrane.

画像



15% SDS-PAGE analysis of ab208317 (3µg).

Observed MW: 18-28KDa (SDS-PAGE under reducing conditions).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors