

Recombinant Human GTF2E2 protein ab126665

画像数 1

製品の詳細

| | |
|-------------|--|
| 製品名 | Recombinant Human GTF2E2 protein |
| 精製度 | > 90 % SDS-PAGE. Purified using conventional chromatography techniques. |
| 発現系 | Escherichia coli |
| アクセッション番号 | <u>P29084</u> |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MGSSHHHHHH SSGLVPRGSH MGSMDPSLL RERELFKKRA LSTPVVEKRS ASSESSSSSS KKKKTKVEHG GSSGSKQNSD HSNGSFNLKA LSGSSGYKFG VLAKIVNYMK TRHQRGDTHP LTLDEILDET QHLDIGLKQK QWLMTEALVN NPKIEVIDGK YAFKPKYNVR DKKALLRLLD QHDQRGLGGI LLEDIEEALP NSQKAVKALG DQILFVNRPD KKKILFFNDK SCQFSVDEEF QKLWRSVTVD SMDEEKIEEY LKRQGISSMQ ESGPKKVAPI QRRKKPASQK KRRFKTHNEH LAGVLKDYS D ITSSK |
| 予測される分子量 | 36 kDa including tags |
| 領域 | 1 to 291 |
| タグ | His tag N-Terminus |

特性

Our **Abpromise guarantee** covers the use of **ab126665** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|-------------------------------|
| アプリケーション | SDS-PAGE Mass Spectrometry |
| 質量分析 | MALDI-TOF |

製品の状態

Liquid

備考

This product was previously labelled as TFIIIE beta

前処理および保存**保存方法および安定性**

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報**機能**

Recruits TFIIH to the initiation complex and stimulates the RNA polymerase II C-terminal domain kinase and DNA-dependent ATPase activities of TFIIH. Both TFIIH and TFIIIE are required for promoter clearance by RNA polymerase.

配列類似性

Belongs to the TFIIIE beta subunit family.

Contains 1 TFIIIE beta DNA-binding domain.

細胞内局在

Nucleus.

画像

ab126665 (3µg) visualised by SDS-PAGE (15%).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors