

### Recombinant human ICK protein ab101774

画像数 4

#### 製品の詳細

製品名	Recombinant human ICK protein
生理活性	The Specific activity of ab101774 was determined to be 4.8 nmol/min/mg.
精製度	> 75 % SDS-PAGE. Purity was determined to be >75% by densitometry. Affinity purified.
発現系	Baculovirus infected insect cells
アクセッション番号	<b><u>Q9UPZ9</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	98 kDa including tags
領域	1 to 632

#### 特性

Our **Abpromise guarantee** covers the use of **ab101774** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot Functional Studies SDS-PAGE
製品の状態	Liquid
備考	<b><u>ab64311</u></b> (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity

#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
------------	--


This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 関連情報


**関連性** ICK may play a key role in the development of multiple organ systems and particularly in cardiac development.

**細胞内局在** Cytoplasmic and Nuclear


## 画像



The specific activity of ICK (ab101774) was determined to be 4.8 nmol/min/mg as per activity assay protocol




SDS PAGE analysis of ab101774



SDS-PAGE analysis of ab101774.

SDS-PAGE - Recombinant human ICK protein  
(ab101774)



The Specific activity of ab101774 was determined to be 4.8 nmol/min/mg.

Functional Studies - Recombinant human ICK  
protein (ab101774)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors