

# Recombinant Human IKZF3 protein ab132618

画像数 1

### 製品の詳細

---

製品名	Recombinant Human IKZF3 protein
発現系	Wheat germ
アクセッション番号	<b>Q9UKT9</b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MEDIQTNAELKSTQEQSVPAESA AVLNDYSLTKSHEMENVDS GEGPAN EDEDIGDDSMKVKDEYSERDENVLKSEPMGNAEEPEIPYSYS REYNEYEN IKLERHVVSFDSSRPTSGKMNC DVCGLSCISFNVLMVHKRSH TGERPFQC NQCGASFTQGNLLRHIKLHTGEKPFKCHLCNYACRRDALT GHLRTHSV EKPYKCEFCGRSYKQRSSLEE HKERCRTFLQSTDPGDTASAE ARHIKAEM GSERALVLDRLASNVAKRKSSMPQKFIGEKRHCFD VNYNSSY MYEKESEL IQTRMMDQAINNAISYLGAELRPLVQTPAPTSEMVPVISS MYPIALTR AEMSNGAPQELEKKS IHLPEKSVPSERGLSPNNSGHDSTDTD SNHEERQN HIYQQNHMVL SRARNGMPLLKEVPRSYELLKPPPICPRDSVK VINKEGEV MDVYRCDHCRVLF LDYVMFTIHMGC HGRDPPFCNMCGYRSH DRYEFSSH IARGEHRALLK
予測される分子量	82 kDa including tags
領域	1 to 509

### 特性

---

Our **Abpromise guarantee** covers the use of **ab132618** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------	--

#### 関連情報

機能	Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.
組織特異性	Expressed most strongly in peripheral blood leukocytes, the spleen, and the thymus.
配列類似性	Belongs to the Ikaros C2H2-type zinc-finger protein family. Contains 6 C2H2-type zinc fingers.
翻訳後修飾	Phosphorylation on tyrosine residues induced by IL2 is required for dissociation from HRAS and nuclear translocation of IKZF3 in T-cells. Phosphorylation on tyrosine residues induced by IL4 is required for dissociation from Bcl-X(L) in T-cells.
細胞内局在	Nucleus. Cytoplasm.

#### 画像



12.5% SDS-PAGE analysis of ab132618 stained with Coomassie Blue.

SDS-PAGE - Recombinant Human IKZF3 protein  
(ab132618)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors