

Recombinant Human ITIH2/SHAP protein ab127154

製品の詳細

製品名	Recombinant Human ITIH2/SHAP protein
精製度	> 90 % SDS-PAGE. Purified via His tag
発現系	Escherichia coli
アクセッション番号	P19823
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	22 kDa
領域	124 to 321
タグ	His-DHFR tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab127154** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	Protein previously labeled as ITIH2.

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
再構成	Reconstitute with water to desired concentration.

関連情報

機能	May act as a carrier of hyaluronan in serum or as a binding protein between hyaluronan and other
----	--

matrix protein, including those on cell surfaces in tissues to regulate the localization, synthesis and degradation of hyaluronan which are essential to cells undergoing biological processes.

組織特異性

Plasma.

配列類似性

Belongs to the ITIH family.

Contains 1 VIT domain.

Contains 1 VWFA domain.

翻訳後修飾

Heavy chains are linked to bikunin via chondroitin 4-sulfate esterified to the alpha-carboxyl of the C-terminal aspartate after propeptide cleavage.

Phosphorylation sites are present in the extracellular medium.

細胞内局在

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors