

Recombinant Human Kallikrein 5 protein ab151839

製品の詳細

製品名	Recombinant Human Kallikrein 5 protein	
精製度	> 95 % SDS-PAGE. ab151839 has purity greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE. 0.2 µM filtered.	
エンドキシン・レベル	< 1.000 Eu/µg	
発現系	HEK 293 cells	
アクセッション番号	<u>Q9Y337</u>	
タンパク質長	Full length protein	
Animal free	No	
由来	Recombinant	
生物種	Human	
配列	VTEHVLANNVSCDHPSNTVPSGSNQDLGAGAGEDARSDDSS SRIINGS CDMHTQPWQAALLLRPNQLYCGAVLVHPQWLLTAAHCRKKVF RVRLGHYS LSPVYESGQQMFQGVKSIHPHGYSHPGHSNDLMLIKLNRRIR PTKDVRPI NVSSHCPASAGTKCLVSGWTTKSPQVHFVKVLQCLNISVLSQ KRCEAYP RQIDDTMFCAGDKAGRDSQCQDSSGGPVVCNGSLQGLVSWGDI PCARPNRP GVTNLCKFTKWIQETIQANSVDHHHHHH	
予測される分子量	31 kDa including tags	
領域	23 to 293	
タグ	His tag C-Terminus	

特性

Our **Abpromise guarantee** covers the use of **ab151839** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	HPLC
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 5.50

Constituents: 0.39% MES, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

関連情報

機能

May be involved in desquamation.

組織特異性

Expressed in skin, breast, brain and testis. Expressed at the stratum granulosum of palmar skin.

配列類似性

Belongs to the peptidase S1 family. Kallikrein subfamily.

Contains 1 peptidase S1 domain.

細胞内局在

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors