

# Recombinant Human KDELC1 protein ab127809

### 製品の詳細

製品名	Recombinant Human KDELC1 protein
精製度	> 95 % Purified via His tag. Purity determined to be >95% by SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	<b><u>Q6UW63</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	20 kDa
領域	326 to 494
タグ	His-DHFR tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab127809** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
再構成	Reconstitute with water to desired concentration.

### 関連情報

配列類似性	Belongs to the KDELC family. Contains 1 filamin repeat.
-------	--

翻訳後修飾

N-glycosylated.

細胞内局在

Endoplasmic reticulum lumen.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors