

Recombinant Human LECT2 protein ab188467

[2 References](#) [画像数 1](#)

製品の詳細

製品名	Recombinant Human LECT2 protein
精製度	> 95 % Densitometry.
エンドキシン・レベル	< 0.100 Eu/μg
発現系	Escherichia coli
アクセッション番号	<u>O14960</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	GPWANICAGKSSNEIRTCDRHGCGQYSAQRSQRPHQGVDILC SAGSTVYA PFTGMIVGQEKPYQNKNAINNGVRISGRGFCVKMFYIKPIKY KGPIKKGE KLGTLPLQKVYPGIQSHVHIENCSSDPTAYL
予測される分子量	16 kDa including tags
領域	19 to 151
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab188467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry ELISA Western blot
質量分析	LC-MS/MS
製品の状態	Lyophilized

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.

Constituents: 0.24% Tris, 0.29% Sodium chloride

再構成

Add 200 µl of deionized water to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely. Filter sterilize your culture media/working solutions containing this non-sterile product before using in cell culture. Reconstituted protein can be stored at 4°C for a week.

関連情報

機能

Has a neutrophil chemotactic activity. Also a positive regulator of chondrocyte proliferation.

組織特異性

Highly expressed in adult and fetal liver and weakly in testis. Not expressed in bone marrow.


配列類似性

Belongs to the LECT2/MIM-1 family.

細胞内局在

Cytoplasm. Secreted.

画像



14% SDS-PAGE analysis of ab188467.

Lane 1: Reduced and boiled sample, 2.5 µg

Lane 2: Non-reduced and non-boiled sample, 2.5 µg

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors