

# Recombinant Human MAGEA4 protein ab99138

画像数 1

### 製品の詳細

製品名	Recombinant Human MAGEA4 protein
精製度	> 90 % SDS-PAGE. ab99138 was purified using conventional chromatography.
発現系	Escherichia coli
アクセッション番号	<b><u>P43358</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p><b>MGSSHHHHHSSGLVPRGSHMSSEQKSQLHCKPEEGVEAQE</b>  <b>EALGLVGAQA</b>  <b>PTTEEQAADVSSSSPLVPGTLEEVPAAESAGPPQSPQGASAL</b>  <b>PTTISFTC</b>  <b>WRQNEGSSSQEEEGPSTSPDAESLFREALSNKVDELAHFL</b>  <b>RKYRAKEL</b>  <b>VTKAEMLERVIKKNYKRCFPVIFGKASESLKMIFGIDVKEVDP</b>  <b>TSNTYTLV</b>  <b>TCLGLSYDGLLGNNQIFPKTGLLIIVLGTIAMEGDSASEEEI</b>  <b>WEELGVMG</b>  <b>VYDGREHTVYGEPRKLLTQDWVQENYLEYRQVPGSNPARYEF</b>  <b>LWGPRALA</b>  <b>ETSYVKVLEHVVRVNARVRIAYPSLREAALLEEEEGV</b></p>
予測される分子量	37 kDa including tags
領域	1 to 317
タグ	His tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab99138** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry

質量分析	MALDI-TOF
製品の状態	Liquid

#### 前処理および保存

---

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine)
------------	--

#### 関連情報


---

機能	Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.
組織特異性	Expressed in many tumors of several types, such as melanoma, head and neck squamous cell carcinoma, lung carcinoma and breast carcinoma, but not in normal tissues except for testes and placenta.
配列類似性	Contains 1 MAGE domain.

---

#### 画像

---



15% SDS-PAGE analysis of 3µg ab99138

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors