

Recombinant Human MRP2 protein ab112272

画像数 1

製品の詳細

製品名	Recombinant Human MRP2 protein
生理活性	useful for Antibody Production and Protein Array
発現系	Wheat germ
アクセッション番号	<u>Q92887</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	LKGYKRPLTLEDVWEVDEEMKTKTLVSKFETHMKRELQKARR ALQRRQEK SSQQNSGARLPGLNKNQSQSDALVLEDVEKKKKKSGTKKDV PKSWLMKA
予測される分子量	37 kDa including tags
領域	214 to 313
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab112272** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot SDS-PAGE
製品の状態	Liquid
備考	This product is useful for Antibody Production and Protein Array.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------	---


pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

Glutathione is reduced

関連情報

機能	Mediates hepatobiliary excretion of numerous organic anions. May function as a cellular cisplatin transporter.
組織特異性	Found on the apical membrane of polarized cells in liver, kidney and intestine. The highest expression is found in liver.
関連疾患	Defects in ABCC2 are the cause of Dubin-Johnson syndrome (DJS) [MM:237500]. DJS is an autosomal recessive disorder characterized by conjugated hyperbilirubinemia, an increase in the urinary excretion of coproporphyrin isomer I, deposition of melanin-like pigment in hepatocytes, and prolonged retention of sulfobromophthalein, but otherwise normal liver function.
配列類似性	Belongs to the ABC transporter superfamily. ABCC family. Conjugate transporter (TC 3.A.1.208) subfamily. Contains 2 ABC transmembrane type-1 domains. Contains 2 ABC transporter domains.
細胞内局在	Membrane.

画像



12.5% SDS-PAGE Stained with Coomassie Blue

SDS-PAGE - Recombinant Human MRP2 protein
(ab112272)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors