

Recombinant human Neurotrophin 3 protein ab9792

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant human Neurotrophin 3 protein
生理活性	Determined by its ability to bind recombinant human TrkB Fc Chimera in a functional ELISA.
精製度	> 98 % SDS-PAGE. Sterile filtered Greater than 98% pure by HPLC analyses. Endotoxin level is less than 0.1 g/ml (1EU/g).
発現系	Escherichia coli
アクセッション番号	P20783
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MYAEHKSHRG EYSVCDSESL WVTDKSSAID IRGHQVTVLG EIKTGNSPVK QYFYETRCKE ARPVKNGCRG IDDKHWNSQC KTSQTYVRAL TSENNKLVGW RWIRIDTSCV CALSRKIGRT
領域	139 to 257

特性

Our **Abpromise guarantee** covers the use of **ab9792** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE Western blot
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--


再構成

For lot specific reconstitution information please contact our Scientific Support Team.

関連情報

機能	Seems to promote the survival of visceral and proprioceptive sensory neurons.
組織特異性	Brain and peripheral tissues.
配列類似性	Belongs to the NGF-beta family.
細胞内局在	Secreted.

画像



Anti-Neurotrophin 3 antibody ([ab65804](#)) at 1 µg/ml + Recombinant human Neurotrophin 3 protein (ab9792) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ([ab97080](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 30 seconds

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors