

# Recombinant Human OXSR1 protein ab125579

画像数 1

### 製品の詳細

製品名	Recombinant Human OXSR1 protein
精製度	> 95 % Densitometry. Purity was determined to be >95% by densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>O95747</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	86 kDa including tags
領域	1 to 527

### 特性

Our **Abpromise guarantee** covers the use of **ab125579** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE

製品の状態	Liquid
-------	--------

### 前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.29% Sodium chloride
------------	---

### 関連情報

機能	Regulates downstream kinases in response to environmental stress. May also have a function in
----	---

	regulating the actin cytoskeleton.
<b>組織特異性</b>	Ubiquitously expressed in all tissue examined.
<b>配列類似性</b>	Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. STE20 subfamily. Contains 1 protein kinase domain.
<b>翻訳後修飾</b>	Phosphorylated upon DNA damage, probably by ATM or ATR.

## 画像



SDS-PAGE analysis of ab125579.

SDS-PAGE - Recombinant Human OXSR1 protein  
(ab125579)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors