

# Recombinant Human Pallidin protein ab124579

画像数 1

### 製品の詳細

---

製品名	Recombinant Human Pallidin protein
精製度	> 90 % SDS-PAGE. ab124579 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b><u>Q9UL45</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MSVPGPSSPD GALTRPPYCL EAGEPTPGLS DTSPDEGLIE DLTIEDKAVE QLAEGLLSHY LPDLQRSKQA LQELTQNQVV LLDLTLEQEIS KFKECHSMLD INALFAEAKH YHAKLVNIRK EMLMLHEKTS KLLKRALKLQ QKRQKEELER EQQREKEFER EKQLTARPAK RM
予測される分子量	22 kDa including tags
領域	1 to 172
タグ	His tag N-Terminus

### 特性

---

Our **Abpromise guarantee** covers the use of **ab124579** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

### 前処理および保存

---

## 保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 関連情報

### 機能

Involved in the development of lysosome-related organelles, such as melanosomes and platelet-dense granules. May play a role in intracellular vesicle trafficking, particularly in the vesicle-docking and fusion process.

### 組織特異性

Widely expressed.


### 翻訳後修飾

Phosphorylated.

### 細胞内局在

Cytoplasm. Endomembrane system. It can exist as a soluble protein as well as a peripheral membrane protein.

## 画像



15% SDS-PAGE showing ab124579 (3µg).

SDS-PAGE - Recombinant Human Pallidin protein  
(ab124579)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors