

Recombinant Human PD1 protein (Tagged) (PE) ab271668

画像数 2

製品の詳細

製品名	Recombinant Human PD1 protein (Tagged) (PE)
精製度	> 90 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	Q15116
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	LDSPDRPWNP PTFSPALLVV TEGDNATFTC SFSNTSESFV LNWYRMSPSN QTDKLAAPFE DRSQPGQDCR FRVTQLPNGR DFHMSVVRAR RNSGTYLCG AISLAPKAQI KESLRAELRV TERRAEVPTA HPSPSRPAG QFQ
予測される分子量	75 kDa including tags
分子量情報	This protein runs at a higher molecular weight by SDS-PAGE due to glycosylation.
領域	25 to 167
タグ	Fc tag C-Terminus
配列の追加情報	Human IgG1 Fc The protein was randomly labeled with RPE on accessible lysine residues.
標識	PE. Ex: 488nm, Em: 575nm

特性

Our **Abpromise guarantee** covers the use of **ab271668** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

前処理および保存

保存方法および安定性

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

関連情報

機能

Possible cell death inducer, in association with other factors.

関連疾患

Genetic variation in PDCD1 is associated with susceptibility to systemic lupus erythematosus type 2 (SLEB2) [MIM:605218]. Systemic lupus erythematosus is a chronic, inflammatory and often febrile multisystemic disorder of connective tissue. It affects principally the skin, joints, kidneys and serosal membranes. It is thought to represent a failure of the regulatory mechanisms of the autoimmune system.

配列類似性

Contains 1 Ig-like V-type (immunoglobulin-like) domain.


発生段階

Induced at programmed cell death.


細胞内局在

Membrane.

画像



Native PAGE Fluorescence of ab271668.



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors