

Recombinant Human PGRP1B protein ab127456

製品の詳細

製品名	Recombinant Human PGRP1B protein
精製度	> 85 % SDS-PAGE. Purified via His tag
発現系	Escherichia coli
アクセッション番号	<u>Q96LB8</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	23 kDa
領域	30 to 328
タグ	His-DHFR tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab127456** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
再構成	Reconstitute with water to desired concentration.

関連情報

機能	Pattern receptor that binds to murein peptidoglycans (PGN) of Gram-positive bacteria. Has bactericidal activity towards Gram-positive bacteria. May kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. Binds also to Gram-negative bacteria, and has bacteriostatic
----	--

	activity towards Gram-negative bacteria. Plays a role in innate immunity.
組織特異性	Detected in skin epidermis, eccrine sweat glands and ducts, mucous cells in the submandibular salivary gland, mucous cells in the throat, ciliary body epithelial cells of the eye, small intestine, colon, stomach and in mature epithelial cells of the tongue (at protein level). High expression in skin and esophagus. Expressed also to a much lesser extent in the tonsils and thymus.
配列類似性	Belongs to the N-acetylmuramoyl-L-alanine amidase 2 family.
翻訳後修飾	N-glycosylated.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors