

Product datasheet

Recombinant Human PI3 Kinase p110 beta protein ab159128

画像数 1

製品の詳細

製品名	Recombinant Human PI3 Kinase p110 beta protein
発現系	Wheat germ
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	EFRRKMRKFSEEEKILSLVGLSWMDWLKQTYPPPEHEPSIPENL EDKLYGGK LIVAVHFENCQDVFSFQVSPNMNPIKVNELAIQKRLTIHGKE DEVSPYDY VLQVSGRVEY
領域	147 to 256
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab159128** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Western blot ELISA
製品の状態	Liquid
備考	

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------	--

機能	Phosphorylates PtdIns, PtdIns4P and PtdIns(4,5)P2 with a preference for PtdIns(4,5)P2.
組織特異性	Expressed ubiquitously.
パスウェイ	Phospholipid metabolism; phosphatidylinositol phosphate biosynthesis.
配列類似性	Belongs to the PI3/P14-kinase family. Contains 1 PI3K/P14K domain.

画像



ab159128 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors