

# Recombinant Human Plakophilin 1 protein ab127241

### 製品の詳細

製品名	Recombinant Human Plakophilin 1 protein
精製度	> 95 % SDS-PAGE. Purified via His tag
発現系	Escherichia coli
アクセッション番号	<b><u>Q13835</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	29 kDa
領域	54 to 317
タグ	His-DHFR tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab127241** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
再構成	Reconstitute with water to desired concentration.

### 関連情報

機能	Seems to play a role in junctional plaques. Contributes to epidermal morphogenesis.
組織特異性	Isoform 2 is widely expressed. Isoform 1 is expressed in stratified squamous, complex, glandular

	duct and bladder epithelia.
関連疾患	Ectodermal dysplasia-skin fragility syndrome
配列類似性	Belongs to the beta-catenin family. Contains 9 ARM repeats.
細胞内局在	Nucleus and Nucleus. Cell junction, desmosome.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors