

Product datasheet

Recombinant Human pNO40 protein ab104634

画像数 1

製品の詳細

製品名	Recombinant Human pNO40 protein
精製度	> 90 % SDS-PAGE. ab104634 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>Q9NP64</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMGSMNSSGRPETMENLPALYT IFQGEVAMVT DYGAFIKIPGCRKQGLVHRTHMSSCRVDPSEIVDVGDKVWV KLIGREMK NDRIKVLSMKVVNQGTGKDLDPNNVIEEQEERRRRSFQDYT GQKITLEA VLNTTCKKCGCKGHFAKDCFMQPGGTKYSLIPDEEEEEKEEK SAEFKPD PTRNPSRKRKKEKKKKKHRDRKSSSDSDSSDESSTGKRARHT SKDSKAAK KKKKKKKHKKKHKE
予測される分子量	30 kDa including tags
領域	1 to 241
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab104634** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid

備考

Previously labelled as ZCCHC17.

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00174% PMSF, 0.0308% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

配列類似性


Contains 1 CCHC-type zinc finger.

Contains 1 S1 motif domain.

細胞内局在

Nucleus, nucleolus.

画像



15% SDS-PAGE showing ab104634 (3 µg).

Note: Molecular weight on SDS-PAGE will appear higher than the predicted molecular weight of approximately 30 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors