

# Recombinant Human PPIC protein ab151358

### 製品の詳細

---

製品名	Recombinant Human PPIC protein
精製度	> 95 % SDS-PAGE. Supplied as a 0.2 µM filtered solution.
エンドトキシン・レベル	< 1.000 Eu/µg
発現系	Escherichia coli
アクセッション番号	<b><u>P45877</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MSDKIIHLTDDSFDTDVLKADGAILVDFWAEWCGPCKMIAPI LDEIADEY QGKLTVAKLNIDQNPGTAPKYGIRGIPTLLLFKNGEVAATKV GALSKGQL KEFLDANLAGSGSGMHMHHHHSSGLVPRGSGMKETAAAKFE RQHMDSPD LGTDDDDKAMAKRGPSVTAKVFFDVRIKDKDVGRIVIGLFGK VVPKTVEN FVALATGEGYGYKGSKFHRVIKDFMIQGGDITTDGDTGGVS IYGETFPD ENFKLKHYGIGWVSMANAGPDTNGSQFFITLTKPTWLDGKHV VFGKVIDG MTVVHSIELQATD
予測される分子量	34 kDa
領域	31 to 182
タグ	His tag N-Terminus
配列の追加情報	Trx tag on N terminal

### 特性

---

Our **Abpromise guarantee** covers the use of **ab151358** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
----------	----------

製品の状態 Liquid

## 前処理および保存

---

**保存方法および安定性** Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.  
pH: 7.40  
Constituents: 0.88% Sodium chloride, 10% Glycerol (glycerin, glycerine), 89% Phosphate Buffer

## 関連情報

---

**機能** PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

**配列類似性** Belongs to the cyclophilin-type PPlase family.  
Contains 1 PPlase cyclophilin-type domain.

**細胞内局在** Cytoplasm.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors