

# Recombinant Human PR3 protein ab153798

## 1 References [画像数 1](#)

### 製品の詳細

製品名	Recombinant Human PR3 protein	
精製度	> 95 % SDS-PAGE. Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE. Supplied as a 0.2 µM filtered solution.	
エンドトキシン・レベル	< 1.000 Eu/µg	
発現系	HEK 293 cells	
アクセッション番号	<b><u>P24158</u></b>	
タンパク質長	Full length protein	
Animal free	No	
由来	Recombinant	
生物種	Human	
配列	AEIVGGHEAQPHSRPYMASLQMRGNPGSHFCGGTLIHPSFVL TAAHCLRD IPQRLVNVVLGAHNVRTQEPTQQHFSVAQVFLNNYDAENKLN DVLLIQLS SPANLSASVATVQLPQQDQVPHGTQCLAMGWGRVGAHDPPA QVLQELNV TVVTFFCRPHNICTFVPRRKAGICFGDSGGPLICDGIIQGID SFVIWGCA TRLFPDFTRVALYVDWIRSTLRRVDHHHHHH	
予測される分子量	26 kDa including tags	
領域	26 to 249	
タグ	His tag C-Terminus	

### 特性

Our **Abpromise guarantee** covers the use of **ab153798** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	HPLC
製品の状態	Liquid

## 前処理および保存

### 保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.16% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

## 関連情報

### 機能


Polymorphonuclear leukocyte serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) and causes emphysema when administered by tracheal insufflation to hamsters.

### 配列類似性

Belongs to the peptidase S1 family. Elastase subfamily.

Contains 1 peptidase S1 domain.

## 画像



SDS-PAGE analysis using ab153798.

Predicted MW 26 kDa with tag.

SDS-PAGE - Recombinant Human PR3 protein  
(ab153798)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors