

Recombinant Human Psoriasin protein ab105577

1 References [画像数 2](#)

製品の詳細

製品名	Recombinant Human Psoriasin protein
精製度	> 90 % SDS-PAGE. ab105577 was purified using conventional chromatography.
発現系	Escherichia coli
アクセッション番号	<u>P31151</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHHSSGLVPRGSHMSNTQAERSIIGMIDMFHKY TRRDDKIDKP SLLTMMKENFPNFLSACDKKGTNYLADVFEKKDKNEDKKIDF SEFLSLLG DIATDYHKQSHGAAPCSGGSQ
予測される分子量	14 kDa including tags
領域	1 to 101
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab105577** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

組織特異性

Fetal ear, skin, and tongue and human cell lines. Highly up-regulated in psoriatic epidermis. Also highly expressed in the urine of bladder squamous cell carcinoma (SCC) bearing patients.


配列類似性

Belongs to the S-101 family.
Contains 2 EF-hand domains.

細胞内局在


Cytoplasm. Secreted. Secreted by a non-classical secretory pathway.

画像



SDS-PAGE - Recombinant Human Psoriasin protein (ab105577)

15% SDS-PAGE analysis of ab105577 (3ug)



Western blot - Recombinant Human Psoriasin protein (ab105577)

Anti-Psoriasin antibody ([ab83534](#)) at 1 µg/ml + Recombinant Human Psoriasin protein (ab105577) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ([ab97080](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 1 minute

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors