

### Recombinant human REP-2 protein ab91088

#### 製品の詳細

製品名	Recombinant human REP-2 protein
生理活性	Activity: 1 pmol of REP-2 will mediate transfer of 2 pmol of geranylgeraniol by RabGGTase onto 1 pmol of Rab7 in 3 minutes.
精製度	> 90 % SDS-PAGE.
発現系	Saccharomyces cerevisiae
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

#### 特性

Our **Abpromise guarantee** covers the use of **ab91088** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid
備考	Previously labelled as CHML.

#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.20 Constituents: 0.077% DTT, 0.595% HEPES, 0.232% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

#### 関連情報

機能	Binds unprenylated Rab proteins, presents it to the catalytic Rab GGTase dimer, and remains
----	---

bound to it after the geranylgeranyl transfer reaction. The component A is thought to be regenerated by transferring its prenylated Rab back to the donor membrane. Less effective than REP-1 in supporting prenylation of Rab3 family.

#### 配列類似性

Belongs to the Rab GDI family.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors