

Recombinant Human Resistin protein ab186979

製品の詳細

製品名	Recombinant Human Resistin protein
精製度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
エンドトキシン・レベル	< 1.000 Eu/μg
発現系	Escherichia coli
アクセッション番号	<u>Q9HD89</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	KTLCSMEEAINERIQEVAGSLIFRAISSIGLECQSVTSRGDL ATCPRGFA VTGCTCGSACGSWDVRAETTCHCQCAGMDWTGARCCRVQPLE HHHHHH
予測される分子量	11 kDa including tags
領域	19 to 108
タグ	His tag C-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab186979** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE HPLC
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -80°C. pH: 3.00
------------	--

Constituent: 0.12% Acetic acid

再構成 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in ddH₂O. Aliquot the reconstituted solution and store at < -20°C.

関連情報

機能	Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.
組織特異性	Expressed only in fatty tissues.
配列類似性	Belongs to the resistin/FIZZ family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors