

# Recombinant Human SDSL protein ab113400

画像数 1

### 製品の詳細

製品名	Recombinant Human SDSL protein	
精製度	> 90 % SDS-PAGE. Conventional chromatography techniques (anion exchange followed by gel filtration with 20mM Tris pH 7.5, 2mM EDTA).	
発現系	Escherichia coli	
アクセッション番号	<u><b>Q96GA7</b></u>	
タンパク質長	Full length protein	
Animal free	No	
由来	Recombinant	
生物種	Human	
配列	<p><b>MGSSHHHHH SSGLVPRGSH MGSHMDGPVA</b>  <b>EHAKQEPFHV VTPLLESWAL SQVAGMPVFL</b>  <b>KCENVQPSGS FKIRGIGHFC QEMAKKGRHLVCSSGGNAG</b>  <b>IAAAYAARKL GIPATIVLPE STSLQVVQRL</b>  <b>QGEGAEVQLT GKVWDEANLR AQELAKRDGW</b>  <b>ENVPPFDHPL IWKGHASLVQ MGSHMDGPVA</b>  <b>EHAKQEPFHV VTPLLESWAL SQVAGMPVFL</b>  <b>KCENVQPSGS FKIRGIGHFC QEMAKKGRH</b>  <b>ELKAVLRTPP GALVLAVGGG GLLAGVVAGL</b>  <b>LEVGWQHVPI IAMETHGAHC FNAAITAGKL</b>  <b>VTLDPDITSVA KSLGAKTVAA RALECMQVCK</b>  <b>IHSEVVEDTE AVSAVQQLLD DERMLVEPAC</b>  <b>GAALAAIYSG LLRRLQAEGC LPPSLTSVVV</b>  <b>IVCGGNNINS RELQALKTHL GQV</b></p>	
予測される分子量	37 kDa including tags	
領域	1 to 329	
タグ	His tag N-Terminus	

### 特性

Our **Abpromise guarantee** covers the use of **ab113400** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry
	SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.04% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	--

#### 関連情報

機能	Has low serine dehydratase and threonine dehydratase activity.
配列類似性	Belongs to the serine/threonine dehydratase family.

#### 画像



15% SDS-PAGE of SDSL protein: (1) MW ladder, (2) protein loaded 3ug

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors