

Recombinant Human SGT1/ECD protein ab211952

製品の詳細

製品名	Recombinant Human SGT1/ECD protein
精製度	> 90 % SDS-PAGE. This protein was expressed in E. coli as inclusion bodies, refolded using unique “temperature shift inclusion body refolding” technology and chromatographically purified.
発現系	Escherichia coli
アクセッション番号	<u>O95905</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

配列	MASMTGGQQMGRGHHHHHHENLYFQGGFELEETMKLATMED TVEYCLFL IPDESRDSDKHKEILQKYIERIITRFAPMLVPYIWQNQPFNL KYPKGGG VPAHMFVTKFGDNIEDWFIVYVIKQITKEFPELVARIEDN DGEFLLIE AADFLPKWLDPENSTNRVFFCHGELCIIPAPRKSAGSWLPT TPPTIPQA LNIITAHSEKILASESIRAAVNRRIRGYPEKIQASLHRAHCF LPAGIVAV LKQRPRLVAAAVQAFYLRDPIDLRACRVFKTFLPETRIMTSV TFTKCLYA QLVQQRFPDRRSGYRLPPSPDPQYRAHELGMKLAHGFEILC SKCSPHFS DCKKSLVTASPLWASFLESLLKNDYFKGLIEGSAQYRERLEM AENYFQLS VDWPSSLAMSPGEEILTLLQTIPFDIEDLKKEAANLPPEDD DQWLDLSP DQLDQLLQEAVGKKESESVSKEEKEQNYDLTEVSESMKAFIS KVSTHKG ELPREPSEAPITFDADSFNLYFDKILGPRPNESDSDLDDDED FECLDSDD DLDFETHEPGEEASLKGTLNLSYMAQMDQELAHTCISKSF TTRNQVEP VSQTTDNNSEEDSGTGESVMAPVDVLDNLVSNILESYSSQA
----	---

予測される分子量	76 kDa including tags
領域	2 to 644
タグ	His-T7 tag N-Terminus
配列の追加情報	Constructed with 31 aa N-terminal T7 / His / TEV cleavage site Tag. NP_009196.

特性

Our **Abpromise guarantee** covers the use of **ab211952** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as SGT1

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituent: 0.32% Tris HCl Contains NaCl, KCl, EDTA, Sucrose, DTT.
------------	--

関連情報

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors