

Recombinant Human SIGLEC15 protein (Tagged) (Biotin) ab271740

画像数 1

製品の詳細

製品名	Recombinant Human SIGLEC15 protein (Tagged) (Biotin)
精製度	>= 90 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	Q6ZMC9
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	F VRTKIDTTEN LLNTEVHSSP AQRWSMQVPP EVSAEAGDAA VLPCTFTHPH RHYDGPLTAI WRAGEPYAGP QVFRCAAARG SELCQTALSL HGRFLLGNP RRNDLSLRVE RLALADDRRY FCRVEFAGDV HDRYESRHGV RLHVTAAPRI VNISVLPSPA HAFRALCTAE GEPPPALAWS GPALGNSLAA VRSPREGHGH LVTAELPALT HDGRYTCTAA NSLGRSEASV YLFRFHGASG AST
予測される分子量	55 kDa
領域	20 to 263
タグ	Avi tag C-Terminus , Fc tag C-Terminus
配列の追加情報	Fc portion of human IgG1. Extracellular domain.
標識	Biotin

特性

Our **Abpromise guarantee** covers the use of **ab271740** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle. Store In the Dark.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

関連情報

機能

Binds sialylated glycoproteins.

組織特異性

Expressed in macrophage and/or dendritic cells of spleen and lymph nodes.

配列類似性

Belongs to the immunoglobulin superfamily. SIGLEC (sialic acid binding Ig-like lectin) family.


Contains 1 Ig-like C2-type (immunoglobulin-like) domain.

Contains 1 Ig-like V-type (immunoglobulin-like) domain.

細胞内局在

Membrane.

画像



SDS-PAGE analysis of 4 µg ab271740.

SDS-PAGE - Recombinant Human SIGLEC15
protein (Tagged) (Biotin) (ab271740)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors