

Recombinant Human Smac/Diablo protein ab87721

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant Human Smac/Diablo protein |
| 精製度 | > 95 % SDS-PAGE. ab87721 is purified by conventional chromatography techniques. |
| 発現系 | Escherichia coli |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MASMTGGQQM GRGSMVPIA QKSEPHSLSS EALMRRVSL VTDSTSTFLS QTTYALIEAI TEYTKAVYTL TSLYRQYTSL LGKMNSEED EVWQVIIGAR AEMTSKHQEY LKLETTWMTA VGLSEMAAEA AYQTGADQAS ITARNHIQLV KLQVEEVHQL SRKAETKLAE AQIEELRQKT QEEGEERAES EQEAYLRED |

特性

Our **Abpromise guarantee** covers the use of **ab87721** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid


前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
pH: 7.50
Constituent: 0.242% Tris

関連情報

| | |
|--------------|--|
| 機能 | Promotes apoptosis by activating caspases in the cytochrome c/Apaf-1/caspase-9 pathway. Acts by opposing the inhibitory activity of inhibitor of apoptosis proteins (IAP). |
| 組織特異性 | Ubiquitously expressed with highest expression in testis. Expression is also high in heart, liver, kidney, spleen, prostate and ovary. Low in brain, lung, thymus and peripheral blood leukocytes. |
| ドメイン | The mature N-terminus mediates interaction with XIAP. |
| 細胞内局在 | Mitochondrion. Released into the cytosol when cells undergo apoptosis. |

画像



15% SDS-PAGE analysis of 3µg ab87721.

SDS-PAGE - Recombinant Human Smac/Diablo protein (ab87721)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors